


# INDICAID™ Early Detection Pregnancy Test Strips




## INTENDED USE

Rapid chromatographic immunoassay | FDA 510(k) cleared (K240643) | Over-the-counter (OTC) use

The INDICAID™ hCG Early Result Pregnancy Test Strip is a rapid chromatographic immunoassay intended for the qualitative detection of human chorionic gonadotropin (hCG) in urine, as an aid in early detection of pregnancy.


## PRODUCT CONFIGURATIONS

|                                                                           |                                                            |
|---------------------------------------------------------------------------|------------------------------------------------------------|
| Contents                                                                  | 25 Individually packaged test strips, instructions for use |
| Specimen type                                                             | Urine                                                      |
| Tests per kit                                                             | 25                                                         |
| UPC                                                                       | 810113420462                                               |
| Retail GTIN                                                               | 00810113420462                                             |
| Case GTIN                                                                 | 10810113420469                                             |
| Kit dimensions (in)                                                       | 2.5" x 2.5" x 4.75"                                        |
| Kit weight                                                                | .15 lb                                                     |
| Case quantity                                                             | 144 kits / 3600 tests                                      |
| Case dimensions                                                           | Provided upon request                                      |
| Case weight                                                               | Provided upon request                                      |
| Pallet quantity                                                           | Provided upon request                                      |
| Pallet dimensions                                                         | Provided upon request                                      |
| Pallet weight                                                             | Provided upon request                                      |
| <b>PERFORMANCE</b>                                                        |                                                            |
| Time to result                                                            | 3 minutes                                                  |
| <b>Analytical Performance</b>                                             |                                                            |
| Sensitivity                                                               | >99% at 10 mIU/mL                                          |
| Validated for early pregnancy detection up to 6 days before missed period |                                                            |
| <b>REGULATORY STATUS</b>                                                  |                                                            |
| FDA clearance                                                             | FDA 510(k) - cleared (K240643)                             |
| OTC use                                                                   | Non-prescription home testing                              |
| HSA/FSA eligible                                                          | Eligible for HSA/FSA reimbursement                         |
| <b>STORAGE CONDITIONS</b>                                                 |                                                            |
| Shelf life                                                                | 36 months                                                  |
| Storage                                                                   | 35–86°F (2–30°C); do not freeze                            |



### Early Detection

Detects pregnancy up to 6 days before the missed period



### Fast answers

Clear results in 3 minutes



### Accessible

Simple, affordable testing with no instrumentation required



### Discreet

At-home testing with no instruments required

### Authorized population

Intended for over-the-counter (OTC) self-testing by lay users

### Setting

Over-the-counter (OTC), non-prescription home use

### Clinical indication

For early detection of pregnancy, as early as six 6 days before the missed period. Results should be confirmed and followed up with a healthcare provider for clinical decision-making.

## LIMITATIONS

- As with any diagnostic test, a confirmed pregnancy diagnosis should be made by a healthcare provider after evaluating all clinical and laboratory findings
- Negative results may occur if the test is taken too early or if urine is too diluted; repeat testing after a few days is recommended if pregnancy is still suspected
- Positive results should be confirmed by a healthcare provider
- Certain fertility treatments or medical conditions involving hCG may cause false results

## ORDERING INFORMATION

Catalog number: IND-DC-HCG-25PPK

657-233-5880 | [USSales@phasesci.com](mailto:USSales@phasesci.com)  
[phasescientificamericas.com](http://phasescientificamericas.com)

This product has been cleared by the U.S. Food and Drug Administration (FDA) under 510(k) for the detection of human chorionic gonadotropin (hCG) in urine, as an aid in the early detection of pregnancy. Results should be interpreted according to the instructions and followed up with a healthcare provider for confirmation and care guidance.

MKT-MC-0061 V 1.1