

# INDICAID™ Urinary Tract Infection Rapid Test

Peace of mind starts here.

Detect signs of a UTI quickly, from the comfort of home. A fast and easy-to-use at-home test to check for urinary tract infection—no appointment required.



## PRODUCT OVERVIEW

### Urinary Tract Infection Rapid Test

The UTI rapid test is a reliable, over-the-counter diagnostic tool for detecting key indicators of urinary tract infections at home. With simple instructions, clear visual results, and no need for specialized equipment, this test empowers users to take control of their health in just 2 minutes.

OTC

|                    |                                    |
|--------------------|------------------------------------|
| Type               | Urinary tract infection test strip |
| Intended use       | Over-the-counter (OTC) / at-home   |
| Sample type        | Urine sample                       |
| Time to result     | 2 minutes                          |
| Pack size          | 1-test kit / 3-test kit            |
| Regulatory status  | FDA 510(k) cleared, CLIA-waived    |
| HSA / FSA approved | Yes                                |



#### Simple

No equipment or training required




#### 510(k) Cleared

Meets regulatory standards for safety and effectiveness.



#### Accessible

Available OTC, HSA/ FSA eligible



#### Fast

Results in just 2 minutes, empowering users to take action quickly

### Star-rated for reliability



4.3/5

Trusted by Amazon buyers nationwide

Amazon customer ratings are based on individual user experiences and do not constitute clinical performance data. Results may vary.



### Kit contents

- 1 or 3 test strips (based on pack size)
- Color interpretation chart
- Instructions for use

## RETAILER INSIGHTS

### Shelf-ready. Demand-driven. Clinically trusted.



- Growing category  
UTIs affect millions annually, especially women, driving demand for at-home tests
- Everyday utility  
UTI symptoms can occur at any time, which supports consistent, year-round sales
- Expands OTC wellness  
Matches retail trends in self-testing and proactive health management
- Broad appeal  
Relevant for women, older adults, caregivers, and people on the go

#### Convenient

Screen from home,  
anytime

#### Actionable

Clear results support timely  
decisions

#### Private

Discreet solution without a  
trip to the doctor



#### Who it's for

- All ages and genders, including children, at risk of UTIs
- Those prone to UTIs due to anatomy or lifestyle
- Individuals with frequent or recurring urinary tract symptoms
- People managing UTIs at home with provider guidance
- Individuals seeking quick health insights
- Caregivers supporting older adults or catheter users

#### Why list our product

Fast-moving SKU with a  
broad target demographic

Drives frequent turn  
with both repeat  
buyers and first-time  
testers alike

Sized and priced for  
strong performance  
across pharmacy shelves  
and wellness kits

No prescription,  
equipment, or training  
required, just fast  
answers from home

Contact us today to place an order: 657-233-5880 | [USsales@phasesci.com](mailto:USsales@phasesci.com)

This test has been cleared by the U.S. Food and Drug Administration (FDA) under 510(k) for the detection of leukocytes and nitrites in urine as an aid in identifying urinary tract infections. It is intended for over-the-counter use by consumers. Results should be interpreted according to the instructions and followed up with a healthcare provider as appropriate.

MKT-MC-0009 V 2.0